


# F319B



## F-319B Diamond Deamabrasion



1. Power Switch
2. Power Energy Increase +
3. Power Energy Decrease -
4. Time Increase +
5. Time Decrease -
6. Start/ Pause



## Microdermoabrasión


Non-invasive procedure that is performed at the level of the corneal layer of the skin and seeks the renewal of the epidermis through mechanical exfoliation, helping to remove dead skin by sanding and sweeping, producing cell regeneration, increasing collagen production, and helping to improve the texture and color of the skin making it more lush and smooth.

It can be performed on the face, neck, décolleté, arms, hands and back, in order to promote the formation of new cells that will replace those eliminated. At the end of the microdermabrasion you may feel a burning sensation and dryness of the skin in the area that was exposed to the procedure, this is normal, you can place the cool hammer to reduce the redness or a refreshing mask. However, depending on the intensity of the treatment, the skin may take on a redder color that should disappear in the following hours.

After turning on the equipment, the suction intensity will be regulated, and the appropriate diamond tip will be used according to the type of skin.

The skin should be clean and dry, and sun exposure should be avoided at least 1 week before. Refrain from using products with retinoids or glycolic acid in the 10 days before and after.

## Punta de Diamante




# F319B



## F-319B Diamond Deamabrasion



1. Power Switch
2. Power Energy Increase +
3. Power Energy Decrease -
4. Time Increase +
5. Time Decrease -
6. Start/ Pause



## Microdermoabrasión

### How to use:

- Connect the tube to the machine and to the microdermabrasion handpiece.
- Place a foam filter before threading the diamond tip. Remove the filter after each treatment.
- Select the tip you wish to use, then screw it onto the handpiece tip until the rubber O-ring is tightened.
- Adjust the suction intensity with the vacuum regulator according to the following table, maintaining communication with the client for safety.
- Place the handpiece with the appropriate diamond tip in contact with the skin. Proceed with a slow and smooth gliding motion. Do not leave the handpiece unmoved for too long on the same skin surface, as this may cause injury.
- Gently stretch the skin with one hand while sliding the diamond tip 2 or 3 times over the same area. Respect this point to avoid redness or injury to the skin.
- After the microdermabrasion treatment, clean the skin of impurities and dead skin with a soft brush, damp sponge or damp towel.
- Use a brush or cotton lightly moistened with alcohol or a disinfectant to clean and disinfect the diamond tip. It is also recommended to sterilize the tips used after each treatment, in an ultraviolet sterilizer or at high temperature.

## INTENSIDAD RECOMENDADA

AREA TRATAMIENTO	PIEL NORMAL	PIEL SENSIBLE
Cara	12-22	8-12
Frente	12-22	10-12
Pecho	15-20	15-18
Nariz	10-15	10-14
Cuello	18-18	10-15
Área de los ojos	8-12	8-10

# F319B



## Ø 22 MM

D100	Corporal. Indicado para acné, rodillas, codos y áreas con estrías
D180	

## Ø 19 MM

D180	Rostro. Indicado para manchas grandes, poros, envejecimiento por exposición solar
D220	

## Ø 15 MM

D100	Cicatrices, líneas finas y diversas imperfecciones
D200	
D250	

## Ø 8 Y 6 MM

D200	Ojos y ambos lados de la nariz. Arrugas, acné, espinillas y puntos negros
D250	

### Precautions:

- Avoid use on irritated or injured skin areas.
- When decreased suction is noticed, replace the filter.
- Do not use granulation (coarse) tips on delicate areas.
- Clean the tips used on each person.
- Avoid bumps and falls of the tips, as this may cause damage and loss of polishing power.
- The microdermabrasion treatment cannot be applied on sensitive areas such as eyelids, and the eyes must be closed during the treatment, it is recommended to cover them with wet cotton pads.
- It is contraindicated in the following cases: skin tumor, moles, birthmarks, angiomas, eyes, herpes and warts.
- Do not attempt to repair the equipment. In case of failure, contact your distributor.

## Contraindicaciones

- Pregnancy and lactation.
- Very sensitive skin.
- Rosacea.
- Seborrheic dermatitis.
- Recent injuries, skin with wounds, infections, psoriasis, dermatitis or unhealthy skin.
- Shaved or depilated skin.
- In skins that have recently undergone deep chemical peeling, burns or surgery.
- Skin sensitivity alterations.
- Poor healing.
- Previous history of keloids.
- Areas with tattoos on the face or permanent makeup.
- Skin conditions.
- Metallic implants.
- Infectious processes or fungus.
- Severe acne.
- Caution when working on tissues with injected fillers, and in the case of clients who have recently been injected with hyaluronic acid, Botox or any other type of product, you should wait a reasonable time to be treated with Hydrafacial.



Yubitzay Avila  
Customer Training  
T. +34 93 824 9020  
y.avila@weelko.com  
S. Yubitzay Avila

Weelko Barcelona S. L.  
Pol. Ind. La Valldan  
Camí de Garreta, 92  
08600 Berga  
Barcelona – Spain  
T. +34 93 824 9020  
info@weelko.com  
www.weelko.com